

## Important Information for the Installation of PurityChrom MCC plus 5.10.014 Update

**Don't install several versions of PurityChrom MCC Plus on one computer!**

**This Update is a full version and will not only update your existing PurityChrom MCC Plus software installation.** If you install the **update** on the same computer as the previous version, the previous version must be **completely removed** from the computer before the new version can be installed. The existing folder must be deleted before.

**Before you delete the folder, make sure that your data have been copied / moved to a safe place.**

**Please copy or move minimum the following data to a safe place**

*PurityChrom.ini (C:\Windows),*

*PurityChrom.cfg, PurityChrom.lic and UserAdministration.cfg (C:\PurityChrom-MCC-Plus)*

*Visualisation Background (C:\PurityChrom-MCC-Plus\Visualisation\Visualisation Backgrounds),*

*Visualisation (C:\PurityChrom-MCC-Plus\Visualisation\Visualisation Files).*

*If the PurityChrom MCC Plus folder includes data, please copy or move this data to a safe place.*

**Uninstall the PurityChrom MCC plus completely and delete the existing folder (C:/PurityChrom-MCC-Plus).**

**The dongle license works independently from the installed PurityChrom MCC plus version number.**

The installation procedure, described below, is for installing the update on the same computer as the previous version. Instead of this, the software can also be installed on a new computer.

### Installation procedure

1. PurityChrom MCC Plus 5.10.014 operates under Microsoft Windows 10, 8.1, 8, 7, Windows Vista, and Windows XP SP3.
2. Store PurityChrom.lic, PurityChrom.cfg, UserAdministration.cfg and visualizations and all PurityChrom MCC Plus data you have already acquired, on a safe place outside of the PurityChrom MCC Plus installation folder.
3. Uninstall "PurityChrom MCC Plus" from *Control Panel – Uninstall programs* (Windows 10 / 8.1 / 8 / 7) or *Control Panel - Add/Remove programs* (Windows XP).
4. Reboot your computer and delete the existing installation folder (C:\PurityChrom-MCC-Plus). Please check before, that you have already copied or moved your existing data.
5. Start the PurityChrom-MCC-Plus Installer. Run the installation setup as administrator. Don't interrupt the installation.

6. If the installation is finished, move the following data:  
**PurityChrom.ini** to C:\Windows  
**PurityChrom.cfg** to C:\PurityChrom-MCC-Plus  
**PurityChrom.lic** to C:\PurityChrom-MCC-Plus  
**UserAdministration.cfg** to C:\PurityChrom-MCC-Plus  
**Visualisation.bmp** to C:\PurityChrom-MCC-Plus\Visualisation\Visualisation Backgrounds  
**Visualisation.vis** to C:\PurityChrom-MCC-Plus\Visualisation\Visualisation Files
7. Restart the PC.